



INSTRUMENTATION DETAILS			
SENSOR	DESCRIPTION	SETTING	ACTION
S1	HIGH DISCH TEMP	2.4HP-220°F	SHUTDOWN
		3.5HP-250°F	
		6.4HP-380°F	
S2	PRESS TRANSDUCER	0-145PSI	OPERATION
S3	BACKUP PRESS SW.	CUT IN-80PSI	BACKUP OPERATION
		CUT OUT-110PSI	
S4	DP SENSOR	+2° C	ALARM
S5	CO SENSOR	10PPM	ALARM

Revised:	Date:	Drawn:	Date:
		GECA	04/03/2012
		This drawing and the information contained herein remain the property of BeaconMedes and may not be used for other than the purpose for which it is loaned without the expressed written permission from BeaconMedes Engineering.	
Description: DWG PNEU DIA 2.4-6.4HP DX RECIP HTM		Part Number: 4107 8531 50	
Scale: NTS		Rev: 00	
Sheet 1 of 1		DO NOT SCALE THIS DOCUMENT	