

BeaconMedæS launches new Oil-less Rotary Vane Vacuum Systems

BeaconMedæS, part of the Atlas Copco Group, introduces the new Oil-less Rotary Vane Vacuum Systems featuring touchscreen controls with an embedded web server and a space-saving design.



The new BeaconMedæS Oil-less Rotary Vane Vacuum Systems provide the ideal medical vacuum solution for smaller medical facilities and surgery centers due to the compact design, quiet operation and low maintenance requirements.

The new vacuum design features state-of-the-art TotalAlert Embedded controls, which allow users to monitor the system remotely from any web browser on the facility's network and to receive electronic notifications of system alarms and service intervals.

“By adding the TotalAlert Embedded control panel to the new Oil-less Rotary Vane Vacuum Systems, a facility can now easily integrate these systems with other TotalAlert source equipment and medical gas alarms,” said Jim Ausman, product manager. “This allows for a complete management of the medical gas network from any computer in the facility.”

The 100% dry technology of the Oil-less Rotary Vane Vacuum pumps, coupled with self-lubricating vanes, means that maintenance requirements are minimal and easy to perform.

BeaconMedæS – Part of the Atlas Copco Group

BeaconMedæS LLC.	Visitors address:	Telephone: + 1 (803) 817-5600
Mailing address	1800 Overview Drive	Telefax: + 1 (803) 817-5750
City, Country	Rock Hill, SC USA	www.beaconmedaes.com

Press Release from the Medical Gas Solutions Division

The simple, compact design with single stage, air-cooled and direct driven vacuum pumps allows for low maintenance and quiet, long lasting operation.

With multiple configurations available meeting a wide range of capacities, the BeaconMedæS Oil-less Rotary Vane Vacuum Systems ensure the right fit for smaller medical vacuum applications.

With the smallest footprint of any of the BeaconMedæS NFPA compliant medical vacuum packages, these systems provide the quiet, vibration-free operation ideal for facilities with limited equipment room space.

For more information, please contact:

Jim Ausman, Product Manager
Phone: +1(803) 817-5624, email: jim.ausman@beaconmedaes.com

BeaconMedæS is a global leader providing complete medical gas solutions, including medical air, medical vacuum, gas distribution equipment and pipeline components. BeaconMedæS has several assembly facilities plus an extensive sales and service presence throughout the world. In July 2006, BeaconMedæS was acquired by Atlas Copco and has since become the global competence center for medical solutions within Atlas Copco. More information on BeaconMedæS can be found at www.beaconmedaes.com

Atlas Copco is a world leading provider of industrial productivity solutions. The products and services range from compressed air and gas equipment, generators, construction and mining equipment, industrial tools and assembly systems, to related aftermarket and rental. In close cooperation with customers and business partners, and with 136 years of experience, Atlas Copco innovates for superior productivity. Headquartered in Stockholm, Sweden, the Group's global reach spans more than 160 markets. In 2008, Atlas Copco had 34 000 employees and revenues of BSEK 74 (BEUR 7.7). Learn more at www.atlascopco.com.

Medical Gas Solutions is a division within Atlas Copco's Compressor Technique business area. It develops, manufactures and markets medical gas equipment, including medical air and vacuum plants, manifolds and pipeline components as well as filtration solutions worldwide. Products are offered under several brands. Production facilities are situated in China, India, Belgium, United Kingdom, Italy and USA. The divisional headquarter is located in Rock Hill, USA.

BeaconMedæS – Part of the Atlas Copco Group

BeaconMedæS LLC.	Visitors address:	Telephone: + 1 (803) 817-5600
Mailing address	1800 Overview Drive	Telefax: + 1 (803) 817-5750
City, Country	Rock Hill, SC USA	www.beaconmedaes.com