



**Series B Premium**  
Front Load DISS





## Introduction

Building on the most common outlet frame in use today, BeaconMedaes is proud to introduce the new Series B Premium Front Load DISS medical gas outlets. This new design features a primary check assembly that is accessed through the front of the outlet, eliminating the need to remove the outlet from the wall and reducing maintenance time. Reducing maintenance time means realizing cost savings over the life of the outlets.

## Series B Family

The Premium Front Load DISS medical gas outlet belongs to the BeaconMedaes Series B family of latch valves. This means the DISS front latch valve is compatible with the full line of Series B wall and console rough-ins, making replacements and upgrades to the new latch valve easy and affordable.

The Premium Front Load DISS also works with the Series B trim plates, including the patented Series B Hookplate. The Series B Hookplate features a retractable arm that takes the load off the outlet, minimizing costly gas leaks. The Hook Plate is treated with Biomaster, an additive that has shown to reduce bacterial growth.\*

\*Disclaimer: Use of this product does not protect user from disease-bearing and food-borne pathogens.



## Reduced Maintenance Time

This new medical gas outlet dramatically reduces the time spent replacing outlet components, resulting in real savings for your facility. Repairing medical gas outlets typically involves removing the front latch valve from the wall/rough-in assembly, before reaching the components that need replacing. With the Front Load DISS from BeaconMedaes, medical gas outlet repair is a quick and easy service. Insert the specialty tool into the front barrel of the latch valve, turn and remove the primary check assembly. It's that easy.



## BeaconMedaes Quality

BeaconMedaes observes the most extensive test procedure for our Series B medical gas outlets. Each latch valve assembly and rough-in are 100% tested for flow, leaks, connector attachment, and indexing. This process is performed on computer-controlled, high-precision test stands, custom built for this unique process. Ensuring the quality that you've come to expect from BeaconMedaes.

The new Premium Front Load DISS outlet is certified to the test standards of UL 1331 for flow rate, leakage, and structural integrity, and meets all requirements of NFPA 99.



---

*Life is in the details.<sup>®</sup>*

