


| INSTRUMENTATION DETAILS | | | |
|-------------------------|-------------------|--------------------|------------------|
| SENSOR | DESCRIPTION | SETTING | ACTION |
| S1 | BACKUP VAC SWT. | 15" Hg ON THE FALL | BACKUP OPERATION |
| S2 | VACUUM TRANSDUCER | VAC-15PSIG | OPERATION |
| S3 | HIGH DISCH TEMP | 15HP - 350°F | SHUTDOWN |



| | | | |
|---|-------|--|---------------------------------------|
| Revised: | Date: | Drawn: | Date: |
| | | JWO | 01/29/14 |
| | | This drawing and the information contained therein remain the property of BeaconMedics and may not be used for other than the purpose for which it is loaned without the expressed written permission from BeaconMedics Engineering. | |
| Description: PNEUMATIC DIAGRAM PX-EXP-HX 15HP CLAW VAC MOD | | Part Number: 4107 8550 22 | Scale: NTS |
| | | HOP131370 Sheet 1 of 1 | Rev: 00 DO NOT SCALE THIS DOCUMENT |