



CONNECTION TABLE
(3 PLACES)

ITEM	DESCRIPTION	SIZE
TP-1	VAC INLET (10HP)	3" NPT
	VAC INLET (15HP)	4" NPT
TP-2	SEAL WATER INLET	1/2" NPT
TP-3	VAC DISCHARGE	4" NPT
TP-4	OVERFLOW DRAIN	1" NPT

PIPING BY OTHERS: - - - - -
FLEX HOSE: NNNNN

BEACONMEDÆS

TITLE:
LL TPX CAMEL TYPE 2 VACUUM SYSTEM
PIPING DIA - GL-35/10HP & GL-60/15HP

DESIGNED BY: T. FELLER	DATE: 10/25/05	This drawing and the information contained therein remain the property of BeaconMedæs and may not be used for other than the purpose for which it is loaned without the expressed written permission from BeaconMedæs Engineering.
CHECKED BY: T. CONNER	DATE: 10/26/05	
REVISED BY:	DATE:	
SCALE: NTS	D	
DRAWING NO. DIA 25 202		REV. A