



CONNECTION TABLE
(2 PLACES)

ITEM	DESCRIPTION	SIZE
TP-1	VAC INLET	4" NPT
TP-2	SEAL WATER INLET	1/2" NPT
TP-3	VAC DISCHARGE	4" NPT
TP-4	OVERFLOW DRAIN	1" NPT
TP-5	CHILLED WATER INLET	1" NPT
TP-6	CHILLED WATER OUTLET	1" NPT

PIPING BY OTHERS: - - - - -
 FLEX HOSE: ~~~~~

BEACONMEDÆS

TITLE:
 LL DPX CAMEL TYPE 2 VACUUM SYSTEM
 PIPING DIA - GL-80/20HP & GL-100/30HP

DRAWN BY: T. FELLER	DATE: 10/25/05	This drawing and the information contained therein remain the property of BeaconMedæS and may not be used for other than the purpose for which it is loaned without the expressed written permission from BeaconMedæS Engineering.
CHECKED BY: T. CONNER	DATE: 10/26/02	
REVISED BY:	DATE:	DRAWING NO.
SCALE: NTS	C	DIA 25 205
		REV. A