


- | ITEM | DESCRIPTION |
|------|--------------------------|
| 1. | SIGHT GLASS |
| 2. | SOLENOID VALVE |
| 3. | VACUUM ISOLATION VALVE |
| 4. | FLOW CONTROL VALVE |
| 5. | CHECK VALVE |
| 6. | VECTRA VACUUM PUMP |
| 7. | RECIRC. CHECK VALVE |
| 8. | ORIFICE UNION |
| 9. | CAMEL WATER RESERVOIR |
| 10. | CHECK VALVE |
| 11. | STRAINER |
| 12. | ANTI-SIPHON VALVE |
| 13. | DIELECTRIC UNION |
| 14. | VACUUM RELIEF VALVE |
| 15. | 120 GALLON RECEIVER TANK |
| 16. | INLET PORT |
| 17. | SEAL WATER PORT |
| 18. | DISCHARGE PORT |



Revised: CCC	Date: 3/11/13	Drawn: CCC	Date: 6/21/13
		This drawing and the information contained therein remain the property of BeaconMedGas and may not be used for other than the purpose for which it is loaned without the expressed written permission from BeaconMedGas Engineering.	
Description: PNEUMATIC DIAGRAM TX VECTRA CAMEL TYPE 1 3-10 HP		Part Number: 4107 8546 29	
Part Number: HOP130608		Rev: 01	
Sheet 1 of 1		DO NOT SCALE THIS DOCUMENT	