



This drawing and the information contained thereon remain the property of BeaconMedics and may not be used for other than the purpose for which it is loaned without the expressed written permission from BeaconMedics Engineering.



INSTRUMENTATION DETAILS			
SENSOR	DESCRIPTION	SETTING	ACTION
S1	BACKUP VAC SWT.	REF. WIRING DIA.	BACKUP OPERATION
S2	VACUUM TRANSDUCER	VAC-15 PSIG	OPERATION

Revised:	Date:	Drawn:	Date:
		JCL	01/12/18
This drawing and the information contained thereon remain the property of BeaconMedics and may not be used for other than the purpose for which it is loaned without the expressed written permission from BeaconMedics Engineering.		Scale: NTS	
Description: PNEUMATIC DIAGRAM DX 1.5-5HP DRY HTM V 80-120G TAE		Part Number: 4107 8554 39	
Sheet 1 of 1		Rev: 00 DO NOT SCALE THIS DOCUMENT	